

High Temperature Constant Temperature Heating Circulator Water Bath Chiller Circulator For Reaction Bath

Item Number: KHB-3



Introduction

Efficient and reliable, KinTek KHB Heating Circulator is perfect for your lab needs. With a max. heating temperature of up to 300°C, it features accurate temperature control and fast heating.

[Learn More](#)

Model	Temperature (°C)	Accuracy (°C)	Pump (L/Min)	Volume (L)	Heating Power (Kw)	Overall Dimensions (mm)
KHB6	Ambient-300	+/-1°C	15	6	1.6	370*380*440
KHB10				10	2.0	370*380*490
KHB15				15	2.5	370*380*490
KHB20				20	3.4	370*380*550
KHB30				30	4.0	480*520*550
KHB50				50	6.5	480*610*620