

# High-Speed Desktop Centrifuge For Blood Sample Processing

Item Number: KT-HPD



## Introduction

Discover the Kintek benchtop high-speed centrifuge for efficient, stable, and precise sample processing in your lab. Ideal for clinical and research use.

[Learn More](#)

| Parameter              | Specification      |
|------------------------|--------------------|
| Display                | LCD Display        |
| Max. Speed             | 16,500 rpm         |
| Speed Control Accuracy | ±10 rpm            |
| Max. RCF               | 18,757 ×g          |
| Max. Capacity          | 12 × 5 mL          |
| Chamber Diameter       | Φ200 mm            |
| Timer Range            | 1s to 99h 59min    |
| Noise Level            | ≤65 dB             |
| Power Supply           | AC 220V±22V 50Hz   |
| Power Consumption      | 420W               |
| Net Weight             | 19 kg              |
| Dimensions (L×W×H)     | 360 × 285 × 258 mm |