

# Vhp Sterilization Equipment Hydrogen Peroxide H2O2 Space Sterilizer

Item Number: KT-C02



## Introduction

A hydrogen peroxide space sterilizer is a device that uses vaporized hydrogen peroxide to decontaminate enclosed spaces. It kills microorganisms by damaging their cellular components and genetic material.

[Learn More](#)

Voltage	AC 220V50HZ
Power	600W
Volume of spray	750ml/60min
Noisy	≤ 45db (A)
Germicide	8%H2O2 solution
Dimension (mm)(L* W * H)	500X500X1130
N.W	55KG