

Graphite Vacuum Furnace IGBT Experimental Graphitization Furnace

Item Number: GF-02



Introduction

IGBT experimental graphitization furnace, a tailored solution for universities and research institutions, with high heating efficiency, user-friendliness, and precise temperature control.

[Learn More](#)

Product model specifications	GF-02-Φ10×15	GF-02-Φ20×30	GF-02-Φ30×40
Volume(L)	1.1	10	28
Limit temperature(C)	3100	3100	3100
Effective heating area (mm)	Φ100×150	Φ200×300	Φ300×400
Power(KW)	30	50	80
Frequency(HZ)	4000	2500	2500
Temperature control method	Japan Shima Electric Thermostat		
Heating method	Induction heating		
Vacuum system	Rotary vane vacuum pump (for high vacuum requirements, Roots vacuum pump and oil diffusion pump are required)		
Sintering atmosphere	N2, Ar		
Rated power supply voltage (V)	380		
Rated heating voltage (V)	According to the design determination, configure the transformer		
Vacuum limit (Pa)	100 (vacuum cold state)		