

Kbr Pellet Press 2T

Item Number: PCKBR



Introduction

Introducing the KINTEK KBR Press - a handheld laboratory hydraulic press designed for entry-level users.

[Learn More](#)

| | |
|---------------------|--------------------|
| Model | PCKBR-2T |
| Working pressure(T) | 0-2(30Mpa) |
| Piston diameter | ≤1MPa/10min |
| Workbench diameter | Φ45mm(D) |
| Number of columns | Two |
| Working space | 54×55 (M×N) |
| Dimensions | 100×220×220(W×L×H) |
| Weight(kg) | 4.8kg |