

Lab Infrared Press Mold

Item Number: PMID



Introduction

Easily release samples from our lab infrared press mold for accurate testing. Ideal for battery, cement, ceramics, and other sample preparation research. Customizable sizes available.

[Learn More](#)

Instrument model	PMID
Sample shape	
Die material	Tungsten carbide
Indenter hardness	HRC68-HRC85
Sample size	Φ13 mm(M)
Cavity depth	20mm(N)
Dimensions	Φ43*78mm(L*H)
Weight	0.76Kg
Diagram of hydraulic powder press size	