

Laboratory Benchtop High Speed Small Centrifuge For Low Temp Serum Separation Dehydration Research

Item Number: KT-HSS



Introduction

Kintek Benchtop High-Speed Centrifuge: Powerful 21,000rpm, precise temperature control. Ideal for clinical, research labs.

[Learn More](#)

Parameter	Specification
Display Mode	LCD Display
Max Speed	21000rpm
Speed Control Accuracy	±20rpm
Max RCF	30642xg
Max Capacity	4 x 750ml
Centrifuge Chamber Specs	Φ420mm
Timer Range	1-99min
Working Noise	≤65dB
Power Supply	AC 220±22V 50Hz
Power	750W
Net Weight	62kg
Product Dimensions (L×W×H)	565×465×395mm