

Laboratory Disc Rotating Mixer

Item Number: KTL-7



Introduction

The laboratory disc rotary mixer can rotate samples smoothly and effectively for mixing, homogenizing and extraction.

[Learn More](#)

Model	LC-TD
Motion mode	Rotation
Voltage	100-240VAC
Frequency	50/60Hz
Power	40W
Angle	0-90°
Motor type	DC motor
Speed range	10-70RPM
Time setting range	1-1199min
Operation mode	Timed/continuous operation
Product size	D280×W210×H300mm
Net weight of product	8Kg
Permissible ambient temperature	5~40°C
Permissible relative humidity	80%RH
Enclosure protection level	IP21