

Laboratory Disc Rotary Mixer For Efficient Sample Mixing And Homogenization

Item Number: KTL-7



Introduction

Efficient Laboratory Disc Rotary Mixer for Precise Sample Mixing, Versatile for Various Applications, DC Motor and Microcomputer Control, Adjustable Speed and Angle.

[Learn More](#)

Model	LC-TD
Motion mode	Rotation
Voltage	100-240VAC
Frequency	50/60Hz
Power	40W
Angle	0-90°
Motor type	DC motor
Speed range	10-70RPM
Time setting range	1-1199min
Operation mode	Timed/continuous operation
Product size	D280×W210×H300mm
Net weight of product	8Kg
Permissible ambient temperature	5□ 40°C
Permissible relative humidity	80%RH
Enclosure protection level	IP21