

Custom Machined And Molded Ptfе Teflon Parts Manufacturer For Laboratory Ito Fto Conductive Glass Cleaning Flower Basket

Item Number: PTFE-13



Introduction

PTFE cleaning racks are mainly made of tetrafluoroethylene. PTFE, known as the "King of Plastics", is a polymer compound made of tetrafluoroethylene.

[Learn More](#)

Diameter*groove width*groove length*applicable capacity	
60mm*1.5mm*16mm*250ml beaker/20 pcs	60mm*2.5mm*16mm*250ml beaker/20 pcs
60mm*1.5mm*21mm*250ml beaker/20 pcs	60mm*2.5mm*21mm*250ml beaker/20 pcs
80mm*1.5mm*26mm*500ml beaker/20 pcs	80mm*2.5mm*26mm*500ml beaker/20 pcs
80mm*1.5mm*31mm*500ml beaker/20 pcs	80mm*2.5mm*31mm*500ml beaker/20 pcs
90mm*2.5mm*41mm*1000ml beaker/22 pcs	