

## Laboratory Ito/Fto Conductive Glass Cleaning Flower Basket

Item Number: PTFE-13



### Introduction

PTFE cleaning racks are mainly made of tetrafluoroethylene. PTFE, known as the "King of Plastics", is a polymer compound made of tetrafluoroethylene.

[Learn More](#)

Diameter\*groove width\*groove length\*applicable capacity

60mm\*1.5mm\*16mm\*250ml beaker/20 pcs

60mm\*2.5mm\*16mm\*250ml beaker/20 pcs

60mm\*1.5mm\*21mm\*250ml beaker/20 pcs

60mm\*2.5mm\*21mm\*250ml beaker/20 pcs

80mm\*1.5mm\*26mm\*500ml beaker/20 pcs

80mm\*2.5mm\*26mm\*500ml beaker/20 pcs

80mm\*1.5mm\*31mm\*500ml beaker/20 pcs

80mm\*2.5mm\*31mm\*500ml beaker/20 pcs

90mm\*2.5mm\*41mm\*1000ml beaker/22 pcs