

# Laboratory Vortex Mixer Orbital Shaker Multifunctional Rotation Oscillation Mixer

Item Number: KTL-4



## Introduction

The inching mixer is small in size, mixes quickly and thoroughly, and the liquid is in a vortex shape, which can mix all the test solutions attached to the tube wall.

[Learn More](#)

Model	KTL-P1	KTL-P2	KTL-P3
Speed range	3500rpm	500~3500rpm	3000rpm
Vibration mode	Press and hold	Press and hold + speed adjustment	Circular
Amplitude	4.5mm	4.5mm	4.5mm
Product size	110*110*75mm	145*160*75mm	110*110*75mm