

Laboratory Multifunctional Small Speed-Adjustable Horizontal Mechanical Shaker For Lab

Item Number: KTL-3



Introduction

The laboratory multifunctional speed-regulating oscillator is a constant-speed experimental equipment specially developed for modern bioengineering production units.

[Learn More](#)

| Model | KTL-3A | KTL-4A | KTL-5A | KTL-6A | KTL-8A |
|----------------------------------|------------------------|------------------------|------------------------|---------------------------------|------------------------|
| Display mode | LED display | LED display | LED display | LED display | LED display |
| Speed range | 40~300rpm | 40~300rpm | 40~300rpm | 40~300rpm | 40~300rpm |
| Timing range | 1~9999min | 1~9999min | 1~9999min | 1~9999min | 1~9999min |
| Standard configuration of pallet | universal spring clamp | universal spring clamp | universal spring clamp | upper clamp, lower spring clamp | universal spring clamp |
| Oscillation model | reciprocating | reciprocating | rotating | rotating | rotating |
| Oscillation amplitude | 20mm | 20mm | 20mm | 20mm | 20mm |
| Pallet working size | 520*300mm | 440*270mm | 440*350mm | 610*460mm (double layer) | 700*620mm |
| Weight | 17.3Kg | 13.5Kg | 20Kg | 65Kg | 65Kg |