

Laboratory Multifunctional Small Speed-Adjustable Horizontal Shaker

Item Number: KTL-3



Introduction

The laboratory multifunctional speed-regulating oscillator is a constant-speed experimental equipment specially developed for modern bioengineering production units.

[Learn More](#)

Model	KTL-3A	KTL-4A	KTL-5A	KTL-6A	KTL-8A
Display mode	LED display	LED display	LED display	LED display	LED display
Speed range	40~300rpm	40~300rpm	40~300rpm	40~300rpm	40~300rpm
Timing range	1~9999min	1~9999min	1~9999min	1~9999min	1~9999min
Standard configuration of pallet	universal spring clamp	universal spring clamp	universal spring clamp	upper clamp, lower spring clamp	universal spring clamp
Oscillation model	reciprocating	reciprocating	rotating	rotating	rotating
Oscillation amplitude	20mm	20mm	20mm	20mm	20mm
Pallet working size	520*300mm	440*270mm	440*350mm	610*460mm (double layer)	700*620mm
Weight	17.3Kg	13.5Kg	20Kg	65Kg	65Kg