

Laboratory Small Constant Temperature Heated Magnetic Stirrer Heater And Stirrer

Item Number: KTL-2



Introduction

The Laboratory Small Constant Temperature Heating Magnetic Stirrer is a versatile tool designed for precise temperature control and efficient mixing in various laboratory applications.

[Learn More](#)

Model	KTL-HA	KTL-HB
Display mode	LCD	LCD
Speed range	100~1600rpm	100~2000rpm
Speed control accuracy	±10rpm	±10rpm
Timing range	1~99h59min	1~99h59min
Temperature control range	RT~300°C	Tray RT~350°C Sample RT~250°C
Temperature control accuracy	Tray ±3°C Solution ±0.5°C	□100°C±0.5°C/ □100°C±1°C (Solution)
Maximum stirring volume (H2O)	5L	20L