

Laboratory Test Sieves And Sieving Machines

Item Number: KT-RCZ



Introduction

Precision lab test sieves & sieving machines for accurate particle analysis. Stainless steel, ISO-compliant, 20µm-125mm range. Request specs now!

[Learn More](#)

Parameter	Specification
Sieve diameter	Φ400, Φ300, Φ200 (mm)
Sieve layers	1-5 (including bottom cover)
Feeding amount (one-time)	≤300g
Vibration amplitude	1-5mm
Excitation force	5KN
Power	0.25KW 220V 50HZ
Speed	1440R/min
Dimensions	600×490×370mm
Sieve weight	60kg
Motion trajectory	circular motion
Timing time	0.01-99.59min
Noise	≤50dB
Material	304 stainless steel
Packing method	wooden box packaging