

## Li-Air Battery Case For Battery Lab Applications

Item Number: BC-04



### Introduction

Lithium air battery (lithium oxygen battery) dedicated battery box. The positive electrode is punched from the inside out, and the inside is smooth.

[Learn More](#)

Dimensions (Dia x Height)	20mm x 3.2mm
Number of openings	17 holes (if 19 holes are required, the hole diameter is 1.2mm)
Material	SS304 Bottom case with Ø12mm x 1.0mm thick mesh disc (Ø1mm holes on it) and top case with PP (Polypropylene) sealing O-ring
Weight	0.10z (2.8 grams)
Application	Excellent for developing Zinc/Lithium-Air battery