

Liquid Nitrogen Cryogenic Grinder Mill Cryomill With Screw Feeder

Item Number: KT-CG02



Introduction

Discover the Liquid Nitrogen Cryogenic Pulverizer with Screw Feeder, perfect for fine material processing. Ideal for plastics, rubber, and more. Boost your lab efficiency now!

[Learn More](#)

| Parameter | Specification |
|--------------------------------|---------------------------------|
| Model | KT-CG02 |
| Host power | 1.8kw |
| Voltage | 220V/50Hz |
| Host speed | 25000r/min |
| Crushing fineness | 40-350 mesh |
| Batch production | 0.3-3kg/h |
| Sample size | 10g |
| Feed size | ≤5mm |
| Seal ring | Low temperature fluorine rubber |
| Discharge mode | Continuous |
| Feeding motor | 40w |
| Feeding speed | 0-10r/min |
| Feeding method | Screw/Manual |
| Liquid nitrogen tank model | YDS-30 |
| Volume of liquid nitrogen tank | 30L(50L optional) |
| Liquid nitrogen pump model | YDB-50J |
| Liquid nitrogen input | Electric pump/Foot pump |
| Liquid nitrogen output control | Scale valve |
| Crushing mode | High speed impact |
| Discharge mode | No sieve wind selection |
| Weight | 46kg |
| Overall dimension | 750450850mm |