

Liquid Nitrogen Cryogenic Grinding Cryomilling Machine With Screw Feeder For Fine Material Processing

Item Number: KT-CG02



Introduction

Discover the Liquid Nitrogen Cryogenic Pulverizer with Screw Feeder, perfect for fine material processing. Ideal for plastics, rubber, and more. Boost your lab efficiency now!

[Learn More](#)

Parameter	Specification
Model	KT-CG02
Host power	1.8kw
Voltage	220V/50Hz
Host speed	25000r/min
Crushing fineness	40-350 mesh
Batch production	0.3-3kg/h
Sample size	10g
Feed size	≤5mm
Seal ring	Low temperature fluorine rubber
Discharge mode	Continuous
Feeding motor	40w
Feeding speed	0-10r/min
Feeding method	Screw/Manual
Liquid nitrogen tank model	YDS-30
Volume of liquid nitrogen tank	30L(50L optional)
Liquid nitrogen pump model	YDB-50J
Liquid nitrogen input	Electric pump/Foot pump
Liquid nitrogen output control	Scale valve
Crushing mode	High speed impact
Discharge mode	No sieve wind selection
Weight	46kg
Overall dimension	750450850mm