

# Metallographic Specimen Mounting Machine For Laboratory Materials And Analysis

Item Number: KT-XQ



## Introduction

Precision metallographic mounting machines for labs—automated, versatile, and efficient. Ideal for sample prep in research and quality control. Contact KINTEK today!

[Learn More](#)

Parameter	Description
Application	Mounting small, difficult to hold or irregular metallographic and petrographic samples
Purpose	Facilitate cutting, grinding, polishing for metallographic, spectral detection or hardness testing
Control System	Program-controlled operation
Modes Available	1. Bakelite powder 2. Electric jade powder 3. User-defined options
Automation	Program-controlled preparation time enables unattended operation
Ideal For	Laboratories in colleges, universities, and industrial/mining enterprises