

Mini Ss High Pressure Autoclave Reactor For Laboratory Use

Item Number: KV-3H



Introduction

Mini SS High Pressure Reactor - Ideal for medicine, chemical, and scientific research industries. Programmed heating temp and stirring speed, up to 22Mpa pressure.

[Learn More](#)

Model	25ml	50ml	100ml	250ml	500ml
Reactor material	304/316L Stainless Steel				
Work temperature	260°C				
Max.work pressure	22Mpa				
Heating method	Modular heating mantle				
Heatint power	600W		800W		
Stirring method	Magnetic stirring/Machanical magnetrton coupling stirring				
Fitting	Inlet/outlet valves, pressre gauge,thermometer				