

No Demolding Lab Infrared Press Mold For Laboratory Applications

Item Number: PMI



Introduction

Effortlessly test your samples with no demolding required using our lab infrared press mold. Enjoy high transmittance and customizable sizes for your convenience.

[Learn More](#)

Instrument model	PMI-A	PMI-B
Sample shape		
Die material	Tungsten carbide	
Indenter hardness	HRC68-HRC85	
Sample size	Φ13 mm(M)	Φ7 mm(M)
Cavity depth	10mm(N)	5mm(N)
Dimensions	Φ76*50*70mm(L*W*H)	Φ76*30*70mm(L*W*H)
Weight	0.76Kg	0.35Kg
Diagram of hydraulic powder press size		