

# Oxygen Probe To Measure Temperature And Active Oxygen Content In Molten Steel

Item Number: KT-T02



## Introduction

Optimize steelmaking with our high-accuracy oxygen probe. Fast, reliable, and essential for precise oxygen and temperature control. Enhance quality and efficiency today.

[Learn More](#)

Parameter	Medium and High Oxygen Probe (OY-TO-Hi)	Low Oxygen (Fixed Aluminum) Probe (OY-TO-Al)
Measurement Range	10-1500ppm	1-20ppm
Reproducibility (Oxygen Potential Deviation)	±3mV	±3mV
Measurement Success Rate	98%	96%
Test Time	4-8s	4-8s
Temperature Measurement	S/R/B model, temperature deviation ±3°C	S/R/B model, temperature deviation ±3°C