

# Rubber Vulcanizer Vulcanizing Machine Plate Vulcanizing Press For Lab

Item Number: RPM-01



## Introduction

The Plate vulcanizing press is a kind of equipment used in the production of rubber products, mainly used for the vulcanization of rubber products. Vulcanization is a key step in rubber processing.

[Learn More](#)

Model No.	RPM-01A	RPM-01B	RPM-01C	RPM-01D
Control Mode	Manual	Electric	PLC control	Vacuum
Capacity	30/50/100ton or customized			
Temperature range	Room temperature~300°C			
Size of working plates	300*300mm or customized			
Work layer	Single layer heating or integrated cooling and heating			
Cooling mode	Water cooling		water cooling or linear cooling	
(W*D*H) Volume	540*400*1250mm	1100*600*1580mm	1100*600*1580mm	1500*600*1580mm
Weight (kg)	210kg	780kg	780kg	900kg
Power supply	3φ AC380V			