


Polygon Press Mold For Lab

Item Number: PMPD



Introduction

Discover precision polygon press molds for sintering. Ideal for pentagon-shaped parts, our molds ensure uniform pressure and stability. Perfect for repeatable, high-quality production.

[Learn More](#)

Instrument model	PMPD
Sample shape	
Die heating	Alloy tool steel :Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size	5*5□10*10□15*15□20*20 mm (M)
Cavity depth	40mm (N)
Dimensions	Φ53*120mm(L*H)
Weight	1.4Kg
Diagram of hydraulic powder press size	