

Portable Autoclave Sterilization Pressure (Digital Display Automatic Type)

Item Number: KT-C05



Introduction

Portable autoclave sterilization pressure is a device that uses pressure saturated steam to quickly and effectively sterilize items.

[Learn More](#)

| | | |
|---------------------------------------|-----------------|----------------|
| Model | KT-18HDD | KT-24HDD |
| Disinfection volume | 18L φ280×260mm | 24L φ280×390mm |
| Fixed work stress | 0.14-0.16Mpa | |
| Rated working temperature | 126°C | |
| Maximum safety pressure | 0.165Mpa | |
| Heat uniformity of sterilization room | ≤±1°C | |
| Power | AC220V.50Hz/2KW | |
| Transport volume (mm) | 420×420×630 | 420×450×750 |
| Gross weight/net weight | 14/12Kg | 15/14Kg |