

## Custom Ptfе Teflon Parts Manufacturer Ptfе Beaker And Lids

Item Number: PTFE-16



### Introduction

The PTFE beaker is a laboratory container that is resistant to acid, alkali, high and low temperatures and is suitable for temperatures ranging from -200°C to +250°C. This beaker has excellent chemical stability and is widely used for heat treatment samples and volume analysis.

[Learn More](#)

Model	inside diameter/mm	height/mm	weight/g
30ml	37	52	25
50ml	44	54	40
100ml	53	54	60
150ml	60	79	89
200ml	70	79	97
250ml	73	88	134
300ml	79	100	167
500ml	95	110	172
1000ml	115	145	364