


## Custom Ptfе Teflon Parts Manufacturer Corrosion-Resistant Cleaning Rack Flower Basket

Item Number: PTFE-20



### Introduction

The PTFE cleaning rack, also known as the PTFE flower basket cleaning flower basket, is a specialized laboratory tool designed for the efficient cleaning of PTFE materials. This cleaning rack ensures thorough and safe cleaning of PTFE items, maintaining their integrity and performance in laboratory settings.

[Learn More](#)

Model	Size	Applicable
Square	60mm*60mm	500ml beaker
Round	Diameter 30mm	75ml beaker
Round	Diameter 45mm	100ml beaker
Round	Diameter 60mm	250ml beaker

We provide custom cleaning flower baskets in all sizes you need