

## Custom Ptfе Teflon Parts Manufacturer For Ptfе Containers

Item Number: PTFE-11



### Introduction

PTFE container is a container with excellent corrosion resistance and chemical inertness.

[Learn More](#)

Round blood outer diameter*inner diameter*inner depth	Inner size of square dish*inner depth
70mm*60mm*15mm	60*60mm *15mm
50mm*40mm*15mm	80*80mm*15mm
60mm*50mm*15mm	90*90mm*15mm
80mm*70mm*15mm	100*100mm*15mm
90mm*80mm*15mm	110*110mm*15mm
100mm*90mm*15mm	120*120mm*15mm
110mm*100mm*15mm	
120mm*110mm*15mm	
130mm*120mm*15mm	
160mm*150mm*15mm	200*200*120MM about 5L PTFE cylinder
92*83*19mm	10L PTFE cylinder
102*93*20mm	15L PTFE cylinder