


Custom Ptfе Teflon Parts Manufacturer For Culture Dish And Evaporation Dish

Item Number: PTFE-26



Introduction

The PTFE culture dish evaporating dish is a versatile laboratory tool known for its chemical resistance and high-temperature stability. PTFE, a fluoropolymer, offers exceptional non-stick properties and durability, making it ideal for various applications in research and industry, including filtration, pyrolysis, and membrane technology.

[Learn More](#)

Model	Bottom outer diameter(mm)	Mouth Inner Diameter(mm)	Total height(mm)	Bottom thickness(mm)	Wall thickness(mm)
Shallow20mm	20	15	10	5	2.5
Shallow30mm	30	25	10	5	2.5
Shallow40mm	40	30	10	5	5
Shallow50mm	50	40	10	5	5
Shallow60mm	60	50	10	5	5
Shallow80mm	80	70	10	5	5
Shallow90mm	90	80	10	5	5
Shallow100mm	100	90	10	5	5
Shallow120mm	120	110	10	5	5
Deep20mm	20	15	20	5	2.5
Deep30mm	30	25	20	5	2.5
Deep40mm	40	30	20	5	5
Deep50mm	50	40	20	5	5
Deep60mm	60	50	20	5	5
Deep80mm	80	70	20	5	5
Deep90mm	90	80	20	5	5
Deep100mm	100	90	20	5	5
Deep120mm	120	110	20	5	5

Model	Outside length(mm)	Inner side length(mm)	Height(mm)	Bottom thickness(mm)	Wall thickness(mm)
Square40mm	40	30	20	5	5
Square50mm	50	40	20	5	5
Square60mm	60	50	20	5	5
Square70mm	70	60	20	5	5

Square80mm	80	70	20	5	5
Square90mm	90	80	20	5	5
Square100mm	100	90	20	5	5
Square110mm	110	100	20	5	5
Square120mm	120	110	20	5	5