

Custom Ptfе Teflon Parts Manufacturer For Microwave Digestion Tank

Item Number: PTFE-37



Introduction

PTFE digestion tanks are renowned for their exceptional chemical resistance, high-temperature stability, and non-stick properties. These tanks are ideal for harsh laboratory environments. Their low coefficient of friction and inert nature prevent chemical interactions, ensuring purity in experimental results.

[Learn More](#)

Model	Diameter(mm)	Height(mm)
30ml	49.5	109
50ml	55	114
100ml	65	124
200ml	74.5	143