


# Custom Ptfе Teflon Parts Manufacturer For Hollow Cleaning Basket And Rack Carrier

Item Number: PTFE-21



## Introduction

The PTFE hollow cleaning flower basket is a specialized laboratory tool designed for efficient and safe cleaning processes. Made from high-quality polytetrafluoroethylene (PTFE), this basket offers exceptional resistance to acids, alkalis, and organic solvents, ensuring durability and reliability in various chemical environments.

[Learn More](#)

Model	size
round	2 inches
round	3 inches
round	4 inches
round	5 inches
round	6 inches
round	8 inches
square	2 inches
square	3 inches
square	4 inches
square	5 inches
square	6 inches
square	8 inches

Model	Diameter	Number of slots	Groove width	Groove length	Overall thickness	Groove depth	Handle length
500ml Beaker Hollow	80mm	24	2.5mm	31mm	8mm	6mm	100mm
500ml Beaker Hollow	80mm	26	2.5mm	26mm	8mm	6mm	100mm
500ml Beaker Hollow	80mm	28	2.5mm	31mm	8mm	6mm	100mm
500ml Beaker Hollow	80mm	30	1.5mm	31mm	8mm	6mm	100mm
500ml Beaker Hollow	80mm	33	1.5mm	26mm	8mm	6mm	100mm
500ml Beaker Hollow	80mm	38	1.5mm	21mm	8mm	6mm	100mm
500ml Beaker Hollow	80mm	40	2.5mm	16mm	8mm	6mm	100mm
500ml Beaker Hollow	80mm	52	1.5mm	16mm	8mm	6mm	100mm
500ml Beaker Hollow	80mm	61	2.5mm	11mm	8mm	6mm	100mm
500ml Beaker Hollow	80mm	78	1.5mm	11mm	8mm	6mm	100mm

250ml Beaker Hollow We provide customized white, off-white, milky white special-shaped flower baskets, including slot size, number of slots and other indicators can be met.

100ml Beaker Hollow We provide customized white, off-white, milky white special-shaped flower baskets, including slot size, number of slots and other indicators can be met.