

Custom Ptfе Teflon Parts Manufacturer For Ptfе Measuring Cylinder 10/50/100MI

Item Number: PTFE-35



Introduction

PTFE measuring cylinder are a rugged alternative to traditional glass cylinders. They are chemically inert over a wide temperature range (up to 260^o C), have excellent corrosion resistance and maintain a low coefficient of friction, ensuring ease of use and cleaning.

[Learn More](#)

Model	Total Height(mm)	Pipe Diameter(mm)	Weight(g)
5ml	83	15	/
10ml	108	17	42
20ml	108	24	64
50ml	118	34	128
100ml	138	40	158
250ml	147	60	374
500ml	160	79	460
1000ml	180	100	699
2000ml	/	/	□