

# Ptfe Reagent Bottle/Wide -Mouth Bottle/Fine Mouth Bottle/Sample Bottle/High Temperature Ptfe

Item Number: PTFE-24



## Introduction

The PTFE Reagent Bottle, also known as PTFE Chemical Bottles or Teflon Reagent Bottle, is a robust alternative to traditional glass reagent bottles. These bottles are highly resistant to both acids and alkalis, featuring a leak-proof screw cap. Ideal for laboratory use, they offer excellent chemical resistance, high temperature capabilities up to 260°C, and superior durability.

[Learn More](#)

Model	Specification	Outer diameter	Height	Inner diameter
PTFE small mouth bottle	25ml	36mm	81mm	21mm
	50ml	49mm	95mm	20mm
	100ml	51mm	98mm	20mm
	150ml	54mm	108mm	20mm
	250ml	70mm	117mm	35mm
	500ml	81mm	147mm	35mm
	1000ml	100mm	188mm	40mm
	2000ml	145mm	217mm	48mm
PTFE wide mouth bottle	50ml	50mm	92mm	35mm
	100ml	52mm	103mm	35mm
	150ml	54mm	106mm	35mm
	250ml	68mm	114mm	40mm
PTFE screw mouth flat bottom reagent bottle	500ml	81mm	149mm	40mm
	1000ml	100mm	192mm	48mm
	2000ml	144mm	220mm	63mm
	5ml	16mm	55mm	*
	10ml	16mm	108mm	*
	20ml	22mm	106mm	*
	30ml	25mm	124mm	*
	50ml	30mm	125mm	*
100ml	39mm	139mm	*	