

## Ptfe Sieve/Ptfe Mesh Sieve/Special For Experiment

Item Number: PTFE-34



### Introduction

PTFE sieve is a specialized test sieve designed for particle analysis in various industries, featuring a non-metallic mesh woven from PTFE (polytetrafluoroethylene) filament. This synthetic mesh is ideal for applications where metal contamination is a concern. PTFE sieves are crucial for maintaining the integrity of samples in sensitive environments, ensuring accurate and reliable results in particle size distribution analysis.

[Learn More](#)

Model	Diameter/cm	Aperture/mm
PTFE Sieve	20	0.15
	20	1
	20	3
	20	5
	20	9.5
	20	10
	20	22.4
PTFE cover	20	5
PTFE Tray	20	5