

Custom Ptfе Teflon Parts Manufacturer F4 Conical Flask Triangular Flask 50 100 250MI

Item Number: PTFE-31



Introduction

The PTFE triangular flask, also known as a Teflon reagent bottle, is a robust, chemical-resistant alternative to traditional glass bottles, suitable for handling both acids and alkalis. These bottles are unbreakable, lightweight, and feature a leak-proof screw cap, making them ideal for laboratory use.

[Learn More](#)

Model	Height including lid(mm)	Bottle bottom diameter(mm)	Inner diameter(mm)
25ml	88	40	27
50ml	100	50	27
100ml	105	59	27
250ml	150	74	37
500ml	174	105	37
1000ml	220	136	37