

Laboratory Sterilizer Lab Autoclave Pulsating Vacuum Desktop Steam Sterilizer

Item Number: KT-B01



Introduction

The pulsating vacuum desktop steam sterilizer is a compact and reliable device used for rapid sterilization of medical, pharmaceutical, and research items.

[Learn More](#)

Model	KT-XD12DV	KT-XD20DV	KT-XD24DV
Sterilization chamber volume	12L ϕ 200×370mm	18L ϕ 250×360mm	24L ϕ 250×470mm
Fixed work stress	0.22Mpa		
Rated working temperature	134°C		
Sterilization temperature selection	105-134°C		
Timing time selection range	0-99min		
Heat uniformity of sterilization room	$\leq \pm 1^\circ\text{C}$		
Power/supply voltage	AC220V.50Hz/1.8KW	AC220V.50Hz/2KW	AC220V.50Hz/2KW
Overall dimensions (mm)	560×470×400	560×470×400	675×470×400
Packing size (mm)	655×560×495	655×560×495	765×560×495
Gross weight/net weight	48/43Kg	52/47Kg	58/53Kg

Model	KT-35DV	KT-50DV
Sterilization chamber volume	35L ϕ 300×500mm	50L ϕ 340×550mm
Fixed work stress	0.22Mpa	
Rated working temperature	134°C	
Sterilization temperature selection	105-134°C	
Timing time selection range	0-99min	
Heat uniformity of sterilization room	$\leq \pm 1^\circ\text{C}$	
Power/supply voltage	AC220V.50Hz/2.7KW	
Overall dimensions (mm)	550×850×500	600×850×560
Packing size (mm)	690×860×700	690×900×700
Gross weight/net weight	113/107Kg	125/110Kg