

# Pulse Vacuum Lifting Sterilizer

Item Number: KT-Q02



## Introduction

The pulse vacuum lifting sterilizer is a state-of-the-art equipment for efficient and precise sterilization. It uses pulsating vacuum technology, customizable cycles, and a user-friendly design for easy operation and safety.

[Learn More](#)

| Model                                     | KT-60SV   | KT-80SV              | KT-100SV  | KT-120SV           |
|---|---|----------------------|---|--------------------|
| Effective volume of sterilization chamber | 60L $\phi$ 396×490mm  | 80L $\phi$ 396×650mm | 100L 504×404×502mm  | 120L 504×404×590mm |
| Height after lifting                      | 1310  | 1590                 | 1475  | 1565               |
| Work pressure                             | 0.22Mpa   |                      |   |                    |
| Operating temperature                     | 134°C   |                      |   |                    |
| Thermal uniformity                        | $\leq \pm 1^\circ\text{C}$                                      |                      |   |                    |
| Timing selection range                    | 0-99min   |                      |   |                    |
| Temperature selection range               | 105-134°C   |                      |   |                    |
| Power/supply voltage                      | Steam generator power: 3KW Drying belt power: 2.4KW/AC220V 50HZ |                      | Steam generator power: 3.6KW Drying belt power: 2.4KW/AC220V 50HZ |                    |
| Dimensions (mm)                           | 686×556×913   | 686×556×1035         | 840×600×970   | 840×600×1060       |
| Transport volume (mm)                     | 789×656×1013  | 786×656×1135         | 940×610×1000  | 940×610×1075       |
| Gross weight/net weight                   | 210Kg/180Kg   | 230Kg/200Kg          | 300Kg/270Kg   | 325Kg/295Kg        |