

# Ring Press Mold For Lab Applications

Item Number: PMO



## Introduction

Ring Press Dies, also known as Circular Pellet Press Die Sets, are integral components in various industrial and laboratory processes.

[Learn More](#)

Instrument model	PMO	
Sample shape		
Die heating	Alloy tool steel :Cr12MoV	
Indenter hardness	HRC60-HRC62	
Sample size	Φ7-3[Φ10-5[Φ20-10 mm M)	Φ30-10[Φ50-20mm (d)
Cavity depth	40mm (N)	45 (d)
Dimensions	Φ53*120mm (L*H)	Φ72*100mm[Φ88*120(D*L)
Weight	1.4Kg	3.5kg[5kg
Diagram of hydraulic powder press size		