

Ring Press Mold For Lab Applications

Item Number: PMO



Introduction

Ring Press Dies, also known as Circular Pellet Press Die Sets, are integral components in various industrial and laboratory processes.

[Learn More](#)

Instrument model	PMQ	
Sample shape		
Die heating	Alloy tool steel :Cr12MoV	
Indenter hardness	HRC60-HRC62	
Sample size	Φ7-3□ Φ10-5□ Φ20-10 mm M)	Φ30-10□ Φ50-20mm (d)
Cavity depth	40mm (N)	45 (d)
Dimensions	Φ53*120mm (L*H)	Φ72*100mm□ Φ88*120(D*L)
Weight	1.4Kg	3.5kg□ 5kg
Diagram of hydraulic powder press size		