


Single Punch Electric Tablet Press Machine Laboratory Powder Tablet Punching Tdp Tablet Press

Item Number: TPM-03



Introduction

The single-punch electric tablet press is a laboratory-scale tablet press suitable for corporate laboratories in pharmaceutical, chemical, food, metallurgical and other industries.

[Learn More](#)

Adjustable punching pressure	50KN
Adjustable tablet diameter	5mm - 25mm
Adjustable tablet thickness	0.5-8mm
Filling depth	20mm
Production capacity	3600 pieces/hour
Motor power	1.5kW
Wooden box packaging Dimensions	580*500*830mm
Weight including outer wooden box	150kg (including export fumigation-free wooden box packaging)