

Slapping Type Sterile Homogenizer Tissue Mashing Homogenizer Dispenser

Item Number: KTL-6



Introduction

The slapping sterile homogenizer can effectively separate the particles contained in and on the surface of solid samples, ensuring that the mixed samples in the sterile bag are fully representative.

[Learn More](#)

Model	KTL-LC400	KTL-LC400M	KTL-LC400GM
Voltage	AC 220V		
Frequency	50Hz		
Total power	200W	220W	500W
Display mode	4.3-inch touch screen		
Slap speed range	3-12 times/s		
Effective volume	3-400ml		
Disinfection function	/	Yes	
Disinfection wavelength	/	253.7nm	
Timing range tapping time	1S-99h59m59m or continuous operation		
Program editing function	8-segment combination		
Net weight	18kg		
Dimensions	460×250×330mm		