

Special Shape Press Mold For Lab

Item Number: PMT



Introduction

Discover high-pressure special shape press molds for diverse applications, from ceramics to automotive parts. Ideal for precise, efficient molding of various shapes and sizes.

[Learn More](#)

Instrument model	PMT
Sample shape	
Die heating	Alloy tool steel :Cr12MoV
Indenter hardness	HRC60-HRC62
Sample size	Φ6□ Φ8□ Φ10□ Φ15□ Φ20mm(M)
Cavity depth	40mm (N)
Dimensions	Φ53*120mm (L*H)
Weight	1.4Kg
Diagram of hydraulic powder press size	