

# Engineering Advanced Fine Alumina Al<sub>2</sub>O<sub>3</sub> Ceramic Crucible For Laboratory Muffle Furnace

Item Number: KM-C01



## Introduction

Alumina ceramic crucibles are used in some materials and metal melting tools, and flat-bottomed crucibles are suitable for melting and processing larger batches of materials with better stability and uniformity.

[Learn More](#)

0.9ml(45x8x6mm)	9ml(60x30x10mm)	30ml(60x40x20mm)	185ml(80x80x40mm)	13ml(100x30x10mm)	300ml(100x100x40mm)	510ml(120x120x50mm)
2ml(30x20x10mm)	10ml(50x20x20mm)	35ml(60x30x30mm)	75ml(90x60x20mm)	25ml(100x30x15mm)	750ml(100x100x100mm)	520ml(140x140x38mm)
3ml(25x15x15mm)	10ml(60x30x14mm)	40ml(50x50x25mm)	165ml(90x60x40mm)	35ml(100x30x20mm)	330ml(115x115x35mm)	425ml(150x150x30mm)
3.5ml(20x20x17mm)	16ml(60x30x15mm)	57ml(75x75x15mm)	200ml(90x60x45mm)	59ml(100x30x30mm)	90ml(120x60x20mm)	750ml(150x150x50mm)
4ml(50x20x10mm)	20ml(60x30x20mm)	68ml(65x65x25mm)	170ml(90x90x30mm)	47ml(100x40x20mm)	190ml(120x60x40mm)	1290ml(150x150x80mm)
5ml(30x20x17mm)	25ml(50x40x20mm)	100ml(65x65x35mm)	300ml(90x90x45mm)	110ml(100x50x30mm)	205ml(120x80x30mm)	800ml(160x160x45mm)
7ml(50x20x15mm)	26ml(50x25x35mm)	110ml(70x70x30mm)	22ml(100x20x20mm)	150ml(100x100x25mm)	260ml(120x80x40mm)	1460ml(180x180x60mm)
3.5ml(20x20x17mm)	10ml(60x30x14mm)	25ml(100x30x20mm)	125ml(98x72x28mm)			
5.5ml(45x22x13mm)	15ml(30x30x30mm)	30ml(80x40x17mm)	1800ml(380x135x60mm)			
9ml(60x25x15mm)	15ml(60x30x15mm)	50ml(90x60x17mm)				
10ml(50x28x14mm)	22ml(100x25x15mm)	84ml(120x60x18mm)				
9ml(60x25x15mm)	50ml(85x30x40mm)	125ml(98x72x28mm)	1050ml(150x150x60mm)	1800ml(380x135x60mm)		
15ml(60x30x15mm)	84ml(120x60x18mm)	460ml(160x120x35mm)	1290ml(150x150x80mm)	4000ml(225x225x100mm)		
15ml(30x30x30mm)	95ml(120x100x15mm)	710ml(155x120x50mm)	1800ml(150x100x150mm)			
22ml(100x25x15mm)	100ml(100x50x30mm)	900ml(120x120x100mm)	4400ml(200x150x200mm)			
0.9ml(45x8x6mm)	59ml(100x30x30mm)	260ml(100x100x35mm)	1880ml(340x220x35mm)	6900ml(330x330x80mm)		
3ml(25x15x15mm)	84ml(65x65x30mm)	425ml(150x150x30mm)	3120ml(300x150x90mm)	14L(330x330x160mm)		
5ml(30x20x17mm)	170ml(170x53x30mm)	445ml(100x100x60mm)	3300ml(200x200x100mm)			
16ml(40x30x20mm)	200ml(115x90x30mm)	1090ml(140x140x70mm)	3760ml(300x200x80mm)			
35ml(70x30x30mm)	230ml(115x65x45mm)	1260ml(160x120x85mm)	4000ml(235x165x130mm)			
One big (50x40x20) and two small (30x20x17)	One large (90x65x30) four small (50x20x20)	One large (100x30x20) four small (20x20x17)	One large (115x90x30) five small (75x20x20)	One big (120x60x20) and two small (50x40x20)		

One big (80x80x40) and two small (60x30x30)	One big (90x65x30) and two small (50x40x20)	One large (100x30x20) and three small (30x20x17)	One large (115x115x35) and three small (100x30x20)	One big (120x60x20) and two small (100x20x20)
One Big (90x60x20) Four small (50x20x20)	One Large (90x65x30) and three small (52x25x23)	One Large (100x40x20) and three small (φ28x17)	One big (120x30x15) five small (20x20x17)	One big (120x60x20) five small (50x20x20)
One large (90x60x17 ) and three small (50x28x1)	One large (100x30x15) and three small (30x20x10)	One Large (100x40x20) and three small (φ30x30)	One large (115x115x35) and three small (100x30x20)	One large (120x60x40) and three small (40x35x25)
One big (90x60x20) and two small (50x40x20)	One big (100x30x20) five small (φ18x26)	One large (100x100x40) nine small (10ml )	One large (120x35x18) seven small (25x15x15)	One large (120x80x30) six small (20ml )
One big (90x60x20) five small (80x10x10)	One large (100x30x20) and four small (φ20x20)	One Large (115x65x45) and four small (50x25x35)	One large (120x60x18 ) and four small (50x28x14)	One large (120x80x40) six small (15ml)

One large (φ58x60) and seven small (φ16x50)	One large (100x40x20) and three small (φ28x17)	One large (120x60x20) two small (100x20x20 semicircle)
One large (φ70x70) and three small (φ28x58)	One large (100x40x20) and three small (φ30x30)	One large (120x60x40) and three small (40x35x25 semicircle)
One large (100x30x20) and four small (φ20x20)	One large (100x100x40) nine small (10ml regular)	One large (120x80x40) six small (15ml)
One big (100x30x20) five small (φ18x26)	One large (115x65x45) four small (50x25x35)	One large (120x80x30) six small (20ml heightened)

2ml(50x10x10mm)	7ml(50x20x15mm)	68ml(65x65x25mm)	110ml(100x50x30mm)	165ml(90x60x40mm)	1000ml(150x110x80)
2.5ml(70x12x7mm)	15ml(75x20x20mm)	75ml(90x60x20mm)	115ml(65x65x35mm)	170ml(90x90x30mm)	2000ml(160x140x120mm)
4ml(50x20x10mm)	34ml(150x20x20mm)	85ml(120x60x18mm)	120ml(90x65x30mm)	185ml(80x80x40mm)	2600ml(200x200x80mm)
4.5ml (80x10x10mm)	57ml(75x75x15mm)	80ml(120x60x18mm)	135ml(95x68x28mm)	185ml(94x68x38mm)	4800ml(250x250x100mm)
7ml(50x25x10mm)	63ml(150x45x20mm)	90ml(120x60x20mm)	160ml(78x28x107)	520ml(140x140x38mm)	

1.5ml with holes (30x10x9.5mm)	60ml(100x43x26mm)
19ml(100x24x17)	114ml(100x56x35mm)
24ml(73x30x22mm)	175ml(100x84x35mm)
36ml(100x27x27mm)	Other specifications can be customized

0.9ml with lid (45x8x6mm)	26ml with lid (50x25x35mm)	80ml with lid (120x60x18mm)	150ml with lid (100x100x25mm)	330ml with lid (115x115x35mm)	520ml with lid (140x140x38mm)
3ml with lid (25x15x15mm)	40ml with lid (50x50x25mm)	84ml with lid (65x65x30mm)	160ml with lid (78x28x107mm)	370ml with lid (200x100x25mm)	750ml with lid (150x150x50mm)
4ml with lid (50x20x10mm)	59ml with lid (100x30x30mm)	100ml with lid (65x65x35mm)	170ml with lid (90x90x30mm)	425ml with lid (150x150x30mm)	800ml with lid (200x100x50mm)
7ml with lid (50x20x15mm)	60ml with lid (105x105x12mm)	100ml with lid (100x50x30mm)	260ml with lid (100x100x35mm)	445ml with lid (100x100x60mm)	900ml with lid (120x120x100mm)
7ml with lid (50x25x10mm)	68ml with lid (65x65x25mm)	110ml with lid (70x70x30mm)	310ml with lid (105x105x40mm)	510ml with lid (120x120x50mm)	
4ml(50x20x10mm)	26ml(50x25x35mm)	195ml(85x80x40mm)	800ml(200x100x50mm)	3100ml(220x220x80mm)	
5.5ml(100x10x10mm)	30ml(60x40x20mm)	200ml(90x60x45mm)	1100ml(305x90x60mm)	298ml(200x100x30mm)	
10ml(75x18x14mm)	35ml(120x30x15mm)	205ml(120x80x30mm)	1480ml(200x200x50mm)		
16ml(60x30x15mm)	150ml(100x100x25mm)	210ml(85x85x40mm)	2400ml(250x155x80mm)		

20ml(52x25x23mm)	190ml(120x60x40mm)	370ml(200x100x25mm)	2830ml(300x180x70mm)
1050ml with lid (150x150x60mm)	1460ml with lid 180x180x60	2600ml with lid 200x200x80	6900ml with lid (330x330x80mm)
1090ml with lid (140x140x70mm)	1480ml with lid 200x200x50	4400ml with cover (200x150x200mm)	11000ml with cover 330x330x120
1290ml with lid (150x150x80mm)	2830ml with cover 300x180x70	4800ml with cover (250x250x100mm)	14L with cover (330x330x160mm)
4ml with lid 20x20x18mm	3.5ml with lid 20x20x17mm		
800ml(130x120x75mm)	1300ml(130x120x115mm)		