


Square Lab Press Mold For Laboratory Applications

Item Number: PMS



Introduction

Create uniform samples easily with Square Lab Press Mold - available in various sizes. Ideal for battery, cement, ceramics, and more. Custom sizes available.

[Learn More](#)

Model	PMS-A	PMS-B	PMS-C	PMS-D	PMS-E	PMS-F	PMS-G
Material	Cr12MoV						
Indenter hardness	HRC60-HRC62						
Sample size	3×3 /4×4 /5×5 /6×6 /8×8 /10×10 mm	12×12 /15×15 /16×16 /18×18 /20×20 mm	22×22 /25×25 /30×30 mm	32×32 /35×35 /40×40 mm	50×50 /60×60 /70×70 mm		81-150 mm(long side) 151-200 mm(long side)
Cavity depth	20 mm	30 mm	40 mm	45 mm	55/60/65 mm		60 mm 60 mm
External dimensions	φ43×93mm	φ53×120mm	φ73×133mm	φ88×150mm	φ98×150mm/φ118×180mm/φ138×180mm		160×140mm 220×160mm
Weight(Kg)	0.65	1.2	2.4	4.8	7.3/11.4/20		25kg-40kg 45kg-80kg

Other sizes can be customized