


# Thin-Layer Spectral Electrolysis Electrochemical Cell

Item Number: ELCST



## Introduction

Discover the benefits of our thin-layer spectral electrolysis cell. Corrosion-resistant, complete specifications, and customizable for your needs.

[Learn More](#)

Parameter	Details
Volume	200 $\mu$ L
Material	All-quartz
Slit Size	2mm (width) $\times$ 10mm (depth)
External Size	12 $\times$ 12 $\times$ 45mm (L $\times$ W $\times$ H)
Internal Size	10 $\times$ 10mm (L $\times$ W)

Component	K003-01 (Platinum Mesh)	K003-02 (Electrode Clip)
Electrode 1	Platinum mesh	Platinum clip
Electrode 2	Platinum wire	Platinum wire
Electrode 3	Standard Ag/AgCl	Mini Ag/AgCl
Included Parts	Quartz cell + PTFE lid + 3 electrodes	Quartz cell + PTFE lid + 3 electrodes