

Laboratory Sterilizer Lab Autoclave Vertical Pressure Steam Sterilizer For Liquid Crystal Display Automatic Type

Item Number: KT-C04



Introduction

Liquid crystal display automatic vertical sterilizer is a safe, reliable and automatic control sterilization equipment, which is composed of heating system, microcomputer control system and overheating and overvoltage protection system.

[Learn More](#)

Model	KT-35HD	KT-50HD	KT-75HD	KT-100HD
Sterilization chamber volume	35L ϕ 318×450mm	50L ϕ 340×550mm	75L ϕ 400×600mm	100L ϕ 440×650mm
Fixed work stress	0.22Mpa			
Rated working temperature	134°C			
Work Pressure	0.25Mpa			
Thermal Uniformity	$\leq \pm 1^\circ\text{C}$			
Timing Selection Range	0-99min/0-99hour59min			
Temperature Selection Range	1-134°C			
Power/supply voltage	1 X 2.5KW/AC220V 50HZ	2 X 2.5KW/AC220V 50HZ	3 X 2.5KW/AC220V 50HZ	
Dimensions (mm)	470×450×1020	520×490×1150	550×530×1180	580×560×1290
Transport Volume (mm)	540×540×1140	580×580×1250	610×610×1280	640×640×1360
Gross weight/net weight	81Kg/63Kg	94Kg/74Kg	113Kg/93Kg	125Kg/105Kg