

Vertical Pressure Steam Sterilizer (Liquid Crystal Display Automatic Type)

Item Number: KT-C04



Introduction

Liquid crystal display automatic vertical sterilizer is a safe, reliable and automatic control sterilization equipment, which is composed of heating system, microcomputer control system and overheating and overvoltage protection system.

[Learn More](#)

| Model | KT-35HD | KT-50HD | KT-75HD | KT-100HD |
|------------------------------|-----------------------|-----------------------|-----------------------|-----------------|
| Sterilization chamber volume | 35L φ318×450mm | 50L φ340×550mm | 75L φ400×600mm | 100L φ440×650mm |
| Fixed work stress | 0.22Mpa | | | |
| Rated working temperature | 134°C | | | |
| Work Pressure | 0.25Mpa | | | |
| Thermal Uniformity | ≤±1°C | | | |
| Timing Selection Range | 0-99min/0-99hour59min | | | |
| Temperature Selection Range | 1-134°C | | | |
| Power/supply voltage | 1 X 2.5KW/AC220V 50HZ | 2 X 2.5KW/AC220V 50HZ | 3 X 2.5KW/AC220V 50HZ | |
| Dimensions (mm) | 470×450×1020 | 520×490×1150 | 550×530×1180 | 580×560×1290 |
| Transport Volume (mm) | 540×540×1140 | 580×580×1250 | 610×610×1280 | 640×640×1360 |
| Gross weight/net weight | 81Kg/63Kg | 94Kg/74Kg | 113Kg/93Kg | 125Kg/105Kg |