

Laboratory Vertical Water Circulating Vacuum Pump For Lab Use

Item Number: KV-1



Introduction

Looking for a reliable water circulating vacuum pump for your lab or small-scale industry? Check out our Vertical Water Circulating Vacuum Pump with five taps and a larger air sucking amount, perfect for evaporation, distillation, and more.

[Learn More](#)

Detailed material

Ejector	check valve	Connecting pipe	three links	pump body	impeller	Air extraction nozzle
Copper	PP	NR	PP	304 stainless steel	304 stainless steel	PP

Model	SHZ-C
Vacuum port	5
Single port flowrate	10L/Min
Vacuum pressure	0.098Mpa
Pump Flowrate	80L/Min
Pump heading	12M
Pump head material	S5304/PTFE
Reservoir volume	50L
Motor power	370W
Dimensions	440*340*830mm
Weight	37Kg