

# Laboratory Wet Three-Dimensional Vibratory Sieve Shaker Machine

Item Number: KT-VM200



## Introduction

The wet three-dimensional vibrating sieving instrument focuses on solving the sieving tasks of dry and wet samples in the laboratory. It is suitable for sieving 20g - 3kg dry, wet or liquid samples.

[Learn More](#)

Maximum load	3 kg
Screen layers	1-8 floors
Maximum sieve tower mass	4 kg
Amplitude	0.2-3 mm
Vibration frequency	3000 times/min
Screening particle size range	20 um-125 mm
The largest available screen diameter	200mm
Maximum sieve tower height	510mm
Additional function	Comfortable visual screen press plate
Screening method	Dry
Electrical description	100-120V/200-240V AC, 50-60Hz, 1500W
Power port	National standard / European standard / American standard / British standard, etc.
Net weight	37 kg
Dimensions (width, thickness and height)	300*405*730mm
Protection level	IP65
Standard	CE/ISO9001